


promea[®]



ISO 9001:2015
QUALITY MANAGEMENT SYSTEM

ISO 13485:2016
MEDICAL DEVICE QUALITY MANAGEMENT

ISO 15197:2013
MEASUREMENT ACCURACY BLOOD GLUCOSE MONITORING SYSTEMS

Dr. Protech-D

Blood Glucose Meter



info@promea.in


www.promea.in



Dr. Protech – D Blood Glucose Meter

Product Intro:

Promea's Dr. Protech D blood glucose meter aims to provide high quality glucose monitoring at a price that allows you to test as often as necessary. The system will give you accurate results in just a few simple steps. It also allows alternative site testing for forearm and palm in addition to fingertip testing.



Features:

- Specially designed, glucose meter aims to provide high quality glucose monitoring at an affordable price.
- No Coding required for the strips.
- Requires only a 0.5 μ L sample to get results in less than 10 seconds.
- Capillary action strips allow a quick and simple application of blood.
- Easy twin button operation.
- Meter memory stores 100 test results.
- Measurement range : 20 – 800 mg/dL, HCT Range : 30 – 55 %.
- Auto shut off : 2 minutes after last action.
- Battery life : 1000 tests.
- Reliable and user friendly.
- Meter Size: 48mmx82mmx18mm
- Display Size: 35 mm x 30 mm
- Weight : 65 gm
- Biosensor Technology – GOD

✉ bdm@promea.in | ☎ 91549 38441

📍 Manufactured by:

Promea Therapeutics Pvt Ltd.

Plot No. M4, Medical Device Park, Sultanpur (V), Patancheru (M),
Sangareddy (Dist), Telangana - 502319. www.promea.in

promea[®]